

new!

BAILEY

Model 25 Deluxe Taping Table

Model 25 Deluxe Taping Table Features:

- Height choices: 30", 36" or 42", no extra charge.
- 1" thick vinyl upholstery in US 393 Black standard, with your choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Full Rear Shelf and End Storage Shelves.
- Stationary 60° Back support.
- Weight capacity 400 lb.
- Ship weight 150 lb.

**NATA
SILVER
CORPORATE
PARTNER**



**Model 25
Deluxe Taping Table**
*[Shown in US 393 Black, with
custom embroidery]*

**See our full catalog of products at
www.baileymfg.com**