

new!

BAILEY

Bariatric Treatment Tables

Plain Bariatric Treatment Table Features:

- 34"W x 28"H x 78"L; optional sizes available.
- 2" thick top, vinyl upholstery in US 387 Camel, standard.
- Choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Standard H-Brace, Optional Plain Shelf available
- ***Six legs and special frame design allows 1000 lb patient capacity!***
- Ships Fully Assembled!
- Ship weight 200 lbs.



Bariatric Treatment Table
[Shown with OPTIONAL Plain Shelf]

Adjustable Back Bariatric Treatment Table Features:

- Same features as Plain Table, with added Adjustable Back
- Adjustable Back features positive locking in 6 to 8 positions.
- Standard H-Brace, Optional Plain Shelf available
- ***Six legs and special frame design allows 1000 lb patient capacity!***
- Ships Fully Assembled!
- Ship weight 200 lbs.



Adjustable Back Bariatric Treatment Table
[Shown in EX 269 Garnet]

See our full catalog of products at
www.baileymfg.com