

**NEW!**

# BAILEY

## Bariatric Treatment Tables

### ***Plain Bariatric Treatment Table Features:***

- 34"W x 28"H x 78"L; optional sizes available.
- 2" thick top, vinyl upholstery in US 387 Camel, standard.
- Choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Standard H-Brace, Optional Plain Shelf available
- ***Six legs and special frame design allows 1000 lb patient capacity!***
- Ships Fully Assembled!
- Ship weight 200 lbs.



**Bariatric Treatment Table**  
*[Shown with OPTIONAL Plain Shelf]*

### ***Adjustable Back Bariatric Treatment Table Features:***

- Same features as Plain Table, with added Adjustable Back
- Adjustable Back features positive locking in 6 to 8 positions.
- Standard H-Brace, Optional Plain Shelf available
- ***Six legs and special frame design allows 1000 lb patient capacity!***
- Ships Fully Assembled!
- Ship weight 200 lbs.



**Adjustable Back Bariatric Treatment Table**  
*[Shown in EX 269 Garnet]*

See our full catalog of products at  
[www.baileymfg.com](http://www.baileymfg.com)