

new!

BAILEY

Model 378 Hand Therapy Table & Desk

Hand Therapy Table Features:

- Durable 48" x 42" top with 24" x 30" therapy area; 30" high
- English Oak top with Natural finish standard.
- Choice of Dark Walnut wood finish, or specified laminate color at additional cost.
- Optional Drawer Locks
- Ships knocked-down for easy assembly.
- Ship weight 200 lb.



*This specially designed table
doubles as a desk.*

*Narrow width and extra cutouts
allow both patient and therapist
unencumbered access.*

*The handy storage area has dual
drawers, two open shelves
and optional locks.*

*Literature rack with removable
dividers also stores flat
therapy boards.*

Hand Therapy Table & Desk
[Shown in Natural Finish]

**See our full catalog of products at
www.baileymfg.com**