

**new!**


# BAILEY

Model 486

## Back Extension Treatment Table

### **Back Extension Table Features:**

- 24"W x 30"H x 78"L; optional sizes available.
- 1" thick top, vinyl upholstery in US 387 Camel, standard.
- Choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Fully adjustable Pelvic Strap; Full Storage Shelf
- Infinitely Adjustable Gas Spring Back with face slot
- Weight capacity 350 lb.
- Ship weight 140 lb.



**Back Extension  
Treatment Table**  
*[Shown in US 417  
Dusty Jade]*

***This specially equipped table is designed with the busy McKenzie Method practitioner in mind.***

***Narrow width and Face Slot eases manual low back therapy procedures.***

***The Pelvic Strap is fully adjustable to stabilize patients of every stature.***

***The generously rounded corners and 1" thick top ensure the comfort of both therapist and patient!***

**See our full catalog of products at  
[www.baileymfg.com](http://www.baileymfg.com)**