


# BAILEY

## Model 6032HD Heavy Lift Box



### *INSIDE DIMENSIONS*



### ***Model 6032HD Features:***

- Heavy Duty Lift Box, **Inside** Dimensions 16"High x 14" Square
- Solid Core 3/4" Hardwood Plywood, Natural Finish
- Weight Capacity 155 LBS., Tare Weight 19 LBS., Total weight 174 LBS., 175 LBS MAX
- NIOSH standard, two high and two low hand slots, and one center weight post, 10" H, uses standard disk type weight plates with 1" DIA center hole, up to 10 KG or 25LB each plate (capacity: 6 each)
- Specially Reinforced stress points, joints and corners for maximum durability
- Furnished fully assembled

***See our full catalog of Occupational Therapy products at  
[www.baileymfg.com](http://www.baileymfg.com)***