


# BAILEY

## Model 4550 Series BariMatic Electric Hi-Lo Upholstered Mat Tables



***BariMatic Electric Mat Table  
[shown in Royal Blue Soft Vinyl]***



***Foot Switch  
(standard)***

### ***Barimatic Mat Series Features:***

- 20" to 28" height adjust, using Foot Switch (standard).
- Adjustable floor levelers.
- Camel color Soft Vinyl top cover standard, with optional Herculite and soft vinyl colors available at no extra charge.
- All steel, heavy duty powder coated frame with crossbrace for stability.
- Lifting Weight capacity 1000 lb.
- Ship weight 270 lb. to 370 lb., depending on size.
- Some Assembly Required.

MODEL	SIZE
4554 BariMatic with Upholstered Top	4' x 7'
4555 Barimatic with Upholstered Top	5' x 7'
4556 Barimatic with Upholstered Top	6' x 8'