

# BAILEY

## Model 28 Portable Sideline Treatment Table The MATT® Folds & Transports For Field Use

### *Mobile Athletic Treatment Table ( The MATT® ) Features:*

- Dimensions: 23" W x 34" H x 66"L; **weight capacity 1000 LBS!**
- Padded Top available in over 15 colors
- Black Powder-Coated Frame locks in place, open and closed, with spring loaded plunger locks
- May be easily moved using 7 inch ball bearing wheels.
- Standard pressed 2" padding; sewn team stripes, or embroidered logo upholstery available.
- End Panel team logos available.
- Ships UPS, (81 LBS) fully assembled

*Invented/trademarked (patent pending) by Pete Miller, veteran Athletic Trainer*



**The MATT@ Table**  
*Just Fold & GO!*



**The MATT@ Table**  
*[Shown in American Beauty Vinyl]*

